

The Flex™ technology breast pump range

Swing Flex™

For quick, efficient expression when you need it



Swing Maxi Flex™

For that extra bottle in half the time



Freestyle Flex™

The flexibility to pump any time, anywhere



Type of pump	Single electric	Double electric	Double electric
Recommended usage	Occasional	Daily	Daily
Easy-to-use and comfortable	•	•	•
PersonalFit Flex™ shield	•	•	•
Adjustable vacuum strength	•	•	•
Vacuum range	-45mmHg to -245mmHg	-45mmHg to -245mmHg	-45mmHg to -245mmHg
2-Phase Expression®	•	•	•
Overflow protection (closed system)	•	•	•
BPA free	•	•	•
Connects to MyMedela app			•
Backlit display showing session timer, battery status, vacuum level, bluetooth connection etc.			•
USB charger			•
Battery operated	Up to 1 hour (battery not included)	Up to 1 hour (battery not included)	Up to 8 double pumping session build-in rechargeable battery
All-in-one solution for storing, cooling and transporting expressed milk			Cooling bag, storage bottles and carry bag
Weight	200g w/o batteries	220g w/o batteries	>300g incl. rechargeable battery
Size	12x13x5.7cm	12x13x5.7cm	14.8x6.6x4.6cm